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#### INSTRUCTIONS FOR USE, STORAGE AND WARRANTY

Dear Customer,  
Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department.  
ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions.  
If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman®. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.  
Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

#### REGULATIONS

**MD** This article is defined as a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prostheses and Orthoses.

#### INDICATIONS

OS6260 Elastic wrist support: Tendonitis, osteoarthritis, mild contusions, minor injuries and instability, inflammatory processes.  
OS6261 Adjustable elastic wrist support: Wrist injury prevention, tendonitis and mild tenosynovitis.  
OS6210 Epicondylitis elbow brace: Epicondylitis (tennis elbow), epitrochleitis (golfer's elbow).  
OS6230 Elastic elbow support: Epicondylitis (tennis elbow), epitrochleitis (golfer's elbow), tendonitis, mild contusions, inflammatory processes, injuries and instability.  
OS6110 Patellar knee band: Patellar tendonitis. Jumper's knee. Osgood-Schlatter. Chondromalacia.  
OS6211 Short elastic knee support / OS6212 Long elastic knee support: Arthritis, arthrosis, traumas and light instabilities of the knee, inflammatory processes.  
OS6801 Elastic thigh support: Prevention of tears of the quadriceps and abductor muscles as well as of the ischiotibial and adductor muscles.  
OS6804 Elastic calf support: Muscle injuries and tibial protection.  
OS6240 Elastic ankle support: Mild ankle instability, joint effusion and inflammation of traumatic origin as well as that caused by previous diseases, postoperative and post-injury irritation (after sprains), tendomyopathies.  
OS6241 Adjustable elastic ankle support: Mild injuries, sprains, dislocations.  
OV02D500 Compression socks: Improvement of venous return, increasing muscle recovery, reducing fatigue during exercise and decreasing recovery periods.  
OS6706 Sport insoles with offloading: Metatarsalgia and talalgia. Pain in the inner longitudinal arch caused by plantar fasciitis or other tendinopathies. Mild to moderate flat foot. Unstable ankle. Runner's knee. Painful or tired feet.  
All Orliman Sport products are indicated for the prevention of injuries in sports and daily activities.

#### FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended.  
If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

#### ⚠ PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the  symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.

Products marked with the  symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.

#### RECOMMENDATIONS-WARNINGS

The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional.

When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

#### RECOMMENDATIONS FOR STORAGE AND WASHING

When not using the product, store it in the original packaging in a dry place at room temperature. Stick the Velcro to each other (if the orthotic device has them), frequently wash by hand with warm water (30° C max.) and mild soap. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

**EN** The manufacturing date is included in the batch number which can be found on the packaging label as **LOT**, in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

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