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INSTRUCTIONS FOR USE, STORAGE AND WARRANTY

Dear Customer,
Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department. ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions.

If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. To use the warranty, go to the place where you bought the product as they must complete the warranty information table on this instructions sheet. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country. Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

REGULATIONS

[MD] This article is defined as a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prostheses and Orthoses.

INDICATIONS

- Post-surgical treatment of the knee.
- Rehabilitation after injury.
- ACL, LCP, LLE, LLI, MENISCI, PATELLAR...

FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended.

If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

When fitting the product, you must adhere to the following instructions:

- 1-Place the patellar ring or the half-moon pad inside the knee brace using the Velcro. If necessary, change the dials articulated by the spiral springs. Place the knee brace on the leg slightly bent. Check its centering in relation to the ball.
- 2-Place you hand in the finger pocket on the back of the thigh. Adjust the desired compression at the level of the thigh then calf by closing the fasteners using Velcro closures.
- 3-Make sure the compression is suitable and that no pressure point remains.
- 4-If necessary, place the "X" strap for reinforcement ligament or half-moon strap for support of the patella according to the indication of your practitioner.

POSITIONING FOR COLLATERAL LIGAMENTS

- 1-Implementation by front place the hollowed out part of the strap at the back of the knee.
- 2-Implement flaps 1 and 2 wrapping knee. Adjust up to a level comfortable and secure to using Velcro.
- 3-End in folding the straps 3 and 4 around the leg and adjust the velcros until whatever the splint is well in place and the optimal confort..

POSITIONING FOR CROSS LIGAMENTS

- 1-Set up by the front. Position the hollowed out part of the strap around the kneecap.
- 2-Implement flaps 1 and 2 wrapping knee. Adjust up to level comfortable and secure to using Velcro.
- 3-End in folding the straps 3 and 4 around the leg and adjust the velcro until whatever the split is well in place and the optimum confort.

HALF-MOON PAD OS SIDE OF PATELLA

- 1-Fix the Velcro lateral strap half-moon on the injured side of the knee.

- 2-Complete the bet in place by stretching simultaneously the ends of the strap around the patella so that the pad exerts a support on the patella

HALF-MOON PAD ON THE PATELLAR TENDON

- 1-Centre and adjust the position on the strap half-moon under the knee cap.

- 2-Lock his position using Velcro.

- 3-Complete the bet in place by stretching the ends of the strap each side of the calf.

PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the  symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.

Products marked with the  symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.

RECOMMENDATIONS-WARNINGS

The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional.


When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

RECOMMENDATIONS FOR STORAGE AND WASHING

When not using the product, store it in the original packaging in a dry place at room temperature. Stick the Velcro to each other (if the orthotic device has them), frequently wash by hand with warm water (30° C max.) and mild soap. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

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 The manufacturing date is included in the batch number which can be found on the packaging label as x, in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

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